



medemagroup



Serial number: \_\_\_\_ - \_\_\_\_\_ - \_\_\_\_\_ - \_\_\_\_\_

Delivery date: \_\_\_\_\_ Year 20\_\_\_\_\_

**This vehicle was supplied by:** \_\_\_\_\_ **Date:** \_\_\_\_\_ / \_\_\_\_\_

---

## Contents

Symbols.....	4
Warning! .....	4
Contagion!.....	4
Introduction.....	5
Declaration of conformity .....	6
Designations .....	7
Joystick parts .....	8
Using the G90T and G90A .....	9
Safety check Joystick .....	10
General care and maintenance .....	11
Insurance.....	12
Preparations / Adjustments prior to use .....	13
Electric seat adjustment .....	14
Seat rotation.....	15
Serial number.....	16
Driving the Mini Crosser M Joy .....	17
Specific driving situations.....	22
Anti-tilt wheels / stabilisers .....	24
Batteries .....	25
Charging.....	25
Battery disposal.....	27
Storage .....	27
Cleaning .....	27
Changing the wheels.....	28
Fuses .....	29
Brakes.....	30
Disengagement.....	31
Transporting by motor vehicle.....	31
Securing to vehicle floor with belts .....	32
Transporting by plane .....	35
Towing .....	35
Dimensioned drawing, Mini Crosser M Joy 4W.....	36
Introduction to seats - Ergo Standard .....	37
Introduction to seats - Mando .....	39
Troubleshooting .....	40
Troubleshooting the Dynamic .....	41
Programming .....	43
Technical data.....	43
International addresses .....	46

## Symbols



Used in the manual to indicate sections describing situations where extra care is required owing to the risk of personal injury.



Used to indicate sections on electromagnetic compatibility (EMC).

---

## Warning!



For safety reasons the vehicle must not be lent to persons who are not completely familiar with it. The vehicle is designed for one person only.



The Mini Crosser M Joy has been designed for users weighing max. 150 kg.

### Joystick



The joystick control box must not be exposed to extremes of temperature or kept in a damp environment for extended periods.



The joystick control box must not be subjected to heavy knocks.



Do not switch off the control box while driving, except in an emergency, as this may damage the electronics.



For cleaning, use a damp cloth with slightly soapy water. Do NOT allow any water or moisture to enter the control box.

---

## Contagion!

### Note!

The tires can sometimes rub off on floor coverings, particularly linoleum. Medema Production assumes no responsibility in case of contagion.

To prevent this, we recommend that you protect delicate floors with some sort of driving surface.

---

## Introduction

Congratulations on your new Mini Crosser M Joy electric mobility scooter.

You have now taken possession of an electric mobility scooter developed for outdoor driving by active users. It is what is called a Class C vehicle in accordance with the European classification of electric scooters.

For optimum enjoyment of this vehicle - and to avoid breakdowns and accidents - we recommend that you read this User Manual carefully. As a new user, you should pay particular attention to the section entitled "Driving the Mini Crosser M Joy".

The Mini Crosser M Joy is designed for safe travel for at least 10 years, up to a max. of 5,000 hours, provided it is serviced and safety-checked every year, corresponding to 500 hours of operation. The service must be carried out by an authorised workshop.



**IMPORTANT!** For safety reasons it is of the utmost importance that service and safety check intervals are complied with, as this minimises the risk of brake failure and short-circuits in the wiring, which could generate heat and cause a fire.

We offer a wide range of accessories for the Mini Crosser M Joy that can make everyday life easier for you. You are always welcome to contact us for further information on special accessories and adaptations.

Medema Production A/S is not responsible for any damage or injuries caused by inappropriate or unsafe use of the Mini Crosser M Joy.

If you have any further questions about the Mini Crosser M Joy or this User Manual, you are always welcome to get in touch. Our contact details are as follows:

**Medema Production A/S**

Tel: +45 7010 2054

Email: [info@minicrosser.com](mailto:info@minicrosser.com)

Internet: [www.minicrosser.dk](http://www.minicrosser.dk)

NB: Errors and omissions excepted. We reserve the right to update this manual as required.

## **CE Declaration of conformity**

Medema Production A/S hereby declares that:

Machine: Mini Crosser  
Use: (Prescribed use in User Manual)  
Model No: M Joy

Complies with the Medical Device Directive 93/42/EEC

The product is made in accordance with the harmonized standard EN 12184 - Electrically powered wheelchairs, scooters and their chargers.

The product is risk analysed in accordance with the harmonized standard DS/EN ISO 14971:2007-04-10 2. edition - Medical devices - Application of risk management to medical devices.




The Mini Crosser can, for a fee, be taken to the nearest dealer for disposal in accordance with current environmental regulations.

---



Manufacturer: Medema Production A/S  
Address: Enggårdvej 7, DK-7400 Herning  
Tel./Fax +45 7010 2054 +45 9716 8582  
Date: 12.04.2010 Signature: \_\_\_\_\_ 

## Designations

The following designations refers to the explanation later in the instruction book.



## Joystick parts

			
Dynamic G90T		Dynamic G90A	
1A	Start/stop button (push button)	1	Start/stop button
1B	Start/stop button (toggle switch)	2	Display
2	Display	3	Select program
3	Gear down/up	4	Display for seat functions (not available)
4	Display for seat functions (not available)	5	Battery indicator
5	Battery indicator	6	Magnetic lock
6	Magnetic lock	7	Horn
7	Horn	8	Joystick
8	Joystick	9	Lights, hazard warning lights and indicators
9	Light, hazard warning lights and indicators	10	Indicator left/lights
10	Select seat function (toggle switch) (not available)	11	Indicator right/hazard warning lights.

### Extra locking system:



Mini Crosser has besides start / stop (I / O button) an additional locking system to ensure unauthorized use of vehicle, in the form of a small magnetic key that comes with a keychain.

Keep magnetic key to the key icon on the joystick, you will hear a signal, and after a few seconds the scooter can not run. Turn the scooter on now (I / O button), a single red lamp is lit. When you want to use the scooter again, hold the magnetic key back to the key symbol, and now you can drive again.

Should you lose the magnetic key, any magnet can be used.



## Using the G90T and G90A

			
Speed selection with G90T		Speed selection with G90A	
1	To start the scooter: Press button 1A or use toggle switch 1B Wait for the battery indicator to stabilise (3-5 seconds)	1	To start the scooter: Press button 1 and wait for the battery indicator to stabilise (3-5 seconds)
2	Select the maximum speed: Keep pressing button 3 or toggle switch 10 to change the numbers between 1, 2, 3, 4 and 5. 1 is the lowest speed and 5 is the highest. (This may vary. Some scooters have a maximum speed of 3).	2	Select the maximum speed: Keep pressing button 3 to change the numbers between 1, 2, 3, 4 and 5 on the display. Select the speed you want. 1 is the lowest speed and 5 is the highest.
Driving		Driving	
1	Start driving by moving the joystick in the direction you want to travel. The speed depends on how far forward you push the joystick. To brake, return the joystick to the starting position in the centre.	1	Start driving by moving the joystick in the direction you want to travel. The speed depends on how far forward you push the joystick. To brake, return the joystick to the starting position in the centre.
2	It is important to check the battery indicator regularly to prevent the scooter stopping because of flat batteries.	2	It is important to check the battery indicator regularly to prevent the scooter stopping because of flat batteries.

### Note!



If you are driving on a sloping surface with poor grip, e.g. gravel or snow, it is important to brake gently to avoid losing control of the scooter. To brake gently: move the joystick to the central position SLOWLY.

## Safety check Joystick

### Daily safety check:

The electronic system has an integrated safety check which runs up to 100 times per minute. To supplement this check, you should carry out the following regular checks.

- Switch off the electronic system (no lights in the display)
- Check if the joystick is bent
- Check if the joystick is damaged in any other way
- Check that it returns to the central position when you release it

If the check reveals any problems, contact a competent service engineer before using the scooter again.

### Weekly safety check:

Parking brake: This test must be carried out on a flat surface with at least one metre of free space around the scooter.

- Start the scooter and slowly move the joystick forward. There is a clicking sound. (The scooter may start to move in this setting).
- Immediately release the joystick and listen for the clicking sound, which should occur within one second.

Repeat in all directions.

- Check that the rubber bellows around the joystick is intact. This is important, as the bellows prevent moisture getting into the electronic system.
- Check that the control box is properly secured.

If the check reveals any problems, contact a competent service engineer before using the scooter again.

### Monthly check

Check the tyre pressure at least once a month. It should be 2.8 bar.

## General care and maintenance

A Mini Crosser M Joy does not require much maintenance. It should be kept in a generally good condition, however. The following should be checked regularly:

- Tyre pressure (if pneumatic tyres are fitted)
- Tyre wear
- Keep the control panel, the charging socket and the electronics box under the seat dry.
- Battery charging



Never wash the Mini Crosser with a high-pressure cleaner or direct water jet! This could damage the Mini Crosser's electronics.

To keep the Mini Crosser in good condition safety-wise, we recommend the following regular checks:

### Daily:

- Test the indicators and driving lights before using the Mini Crosser in the dark or poor visibility.

### Every three months:

Test the brakes and motor disengagement

With the disengagement lever up, it must not be possible to push the Mini Crosser.

Test the brake disengagement function

When the brake disengagement lever is down, the battery indicator should flash to show an error if the Mini Crosser is turned on. In this case the Mini Crosser must not be able to move when the joystick is activated.

Test the handbrake.

Apply the handbrake for a couple of seconds at low speed. This will ensure that the lever arm and brake shoes do not seize up.

Lubricate the lever arm on the brake hub with acid-free oil - left-hand rear wheel.

### Annual

The Mini Crosser M Joy is designed for safe travel for at least 10 years, up to a max. of 5,000 hours, provided it is serviced and safety-checked every year, corresponding to 500 hours of operation. The service must be carried out by an authorised workshop.



**IMPORTANT!** For safety reasons it is of the utmost importance that service and safety check intervals are complied with, as this minimises the risk of brake failure and short-circuits in the wiring, which could generate heat and cause a fire.

(More information in the service manual).

---

## **Insurance**

In the eyes of the law a Mini Crosser M Joy with a maximum speed of 10 km/h is a cycle, so separate insurance is not required.

Most contents/home insurance policies include third-party liability insurance for cyclists and so also cover Mini Crosser M Joy users.


We recommend that you talk to your insurance company about this when the vehicle is delivered. If necessary, comprehensive insurance will have to be taken out separately.

## Preparations / Adjustments prior to use


### Adjusting the height of the seat

Lift the seat off and remove the rear cover. Press the release button out and adjust the seat tube to the desired position. See the picture on the next page. Check that the seat height is correct and press the release button in again. Move the seat tube up/down until the release clicks into a hole.


Make sure you keep your back straight when lifting the seat, which is very heavy. (See below.)




*Lift the seat off the seat tube.*



*The release button pressed in. Press the button out to release the seat tube.*



*Release button pressed out. Now the seat tube can be moved up/down to the height you wish. Press the release button again and make sure it is fitted in one of the 5 holes in the seat tube.*



Seat tube

Release

*Seat tube in top position. There are five holes for adjusting the seat tube. The release must be in a hole for the seat to be stable.*

---

## Electric seat adjustment

The Mini Crosser can be fitted with electric seat adjustment as an optional extra.

To raise the seat, press the switch up. To lower the seat, press the switch down. If the switch is released, the seat will stop automatically. (See figure below.)



Switch for electric seat adjustment

*Electric seat adjustment. If the switch is pressed up, the seat is raised. If the switch is pressed down, the seat is lowered.*

## Seat rotation

Pull the release lever back. The seat can be rotated 90° to either side. When the lever is released, it engages with the seat and holds it in place at 45° intervals.

Other seats that can be supplied for the Mini Crosser work on similar principles. The release lever is normally mounted on the right, but can be put on the left if so wished.


### IMPORTANT!



The Mini Crosser is most stable when the seat is in its lowest position. Always drive carefully when the seat is raised. Never use the seat adjustment when driving on an uneven surface or in hilly terrain.



Take extra care when lowering the seat on a Mini Crosser using electric seat adjustment. Make sure that nothing is trapped in the space between the seat and the chassis.



*Avoid trapping anything when lowering the seat.*






*Release lever for seat rotation. Release lever for seat forward/back.*

## Serial number

All vehicles have a serial number plate showing the year of production, and serial number. The same serial number can also be found on the front of the User Manual.

Please quote this number when making inquiries about servicing, spare parts, etc.

<b>Mini Crosser®</b>		Enggaardvej 7, DK-7400 Herning	
		TEL: +45 7010 1755	
		E-mail: info@minicrosser.com	
 Varenr. / Art.nr. 5707299285358 Serienummer 12589-09-09		 Max. person vægt / vikt - kg 175	
Serienr.: Produktion / Tilverkning:	12589-09-09 2009	Type: Egen vægt/vikt: Max. km/t:	MJ 152 10



Position of serial number plate on vehicle.



## Driving the Mini Crosser M Joy

### Getting in and out

It is important to learn a safe technique for getting in and out of the Mini Crosser.

In general the following method should be followed:

- Make sure that the Mini Crosser M Joy is off (switch turned to 0)
- Make sure that the brake is on. (Lever for disengaging motor in top position.)
- If necessary, turn the seat through 45° or 90° and make sure that it is locked in position (clicked into place).
- If necessary, raise the armrest.

For some users the assistance of an attendant can be recommended. The attendant should:

- Take care not to injure him/herself when lifting/lowering/supporting the user.
- Make sure that the Mini Crosser is stable and unable to move. Turn off the Mini Crosser and check that the brake is on and the seat has been rotated until it clicks into place at either 45° or 90°.
- Make sure that the seat the user is being moved to is stable.



*Turn the Mini Crosser off, rotate the seat and raise the armrest.*

## Braking systems

There are three braking systems on the Mini Crosser:

- Motor brake - adjusts vehicle speed also when going downhill.
- Magnetic brake - the magnetic brake is automatic and engages when the Mini Crosser stops. In an emergency, the Mini Crosser can be stopped instantaneously by turning of the joystick. Please note that this will cause very sharp braking. The rear wheels will lock.  
**Must NOT be used in the normal course of driving.**  
The brake must never be disengaged mechanically using the disengagement lever on a slope. This function is only designed for use when pushing the Mini Crosser on a flat road.
- Handbrake - intended as an emergency brake and parking brake. It must be operated with caution when driving in slippery conditions and downhill.



*When using the handbrake as a parking brake, lock it in braking position by pressing the button in while applying the brake.  
To release the brake, press the button again.*

*Handbrake lock*

## General safety advice

Make sure that the backrest is upright and the seat is as low as possible.

Positioning belts are recommended if the user is unable to maintain a good driving posture independently.

Adjust your driving to road conditions. Take account of light, traffic and weather. Be particularly careful when driving in the dark or in bad weather, such as rain or snow. Avoid driving on gradients with poor surfaces, such as: snow, ice, new-mown grass, wet grass and wet leaves.



Never drive when under the influence. This applies not only to alcohol, but also to drugs and medicines.

Reduce speed immediately if you feel you are losing control.

Always use the indicators when changing direction.

Check that lights and indicators are working before driving off.  
Use your lights when driving after lighting-up time.



ALWAYS switch the scooter off when it is not in use. I/O button.

## Driving

Even though the Mini Crosser is very stable, it can tip over. Avoid sudden changes of speed and direction when travelling at high speed, on poor surfaces and, not least, on slopes.

For short distances the Mini Crosser can drive up steeper gradients than it has been tested as dynamically stable for. The same applies to driving down such gradients. In such cases there is an increased risk of the Mini Crosser tilting and even tipping over. So be extra careful in following the driving tips given below.

Anti-tilt wheels are recommended for driving in very hilly terrain. (Optional extra)

New users are urged to practise the following in an area where there is no other traffic:

- Set the Mini Crosser to low speed. Drive forwards and backwards. Gradually turn the speed selector up and feel the change in the speed of the Mini Crosser.
- Practise starting and braking. Get used to the Mini Crosser's response time.
- Practise driving in a narrow space, similar to inside a shop or through a door.
- Practise turning, and get a sense of how much space is required. Always drive slowly when turning. Practise reversing too.
- Practise cornering and driving over obstacles and on slopes.

Always drive straight up/down kerbs and ramps. Never at an angle. See the illustrations on the following pages.

- Try braking at different speeds and notice the stopping distances.
- If possible, try driving on different surfaces (road, grass and gravel).
- Practise assessing how far you can drive on a single battery charge. Note how quickly the battery indicator changes from green to amber to red.

**Note!**



The driving distance of the scooter will be reduced when driving in hilly districts, into a head wind, in cold weather and with low tyre pressure.

**Traffic regulations**

The traffic legislation for scooters varies from country to country. Before starting to use the vehicle outdoors, it is the user's responsibility to familiarise him/herself with the relevant legislation.

**The following rules apply in Denmark:**

Mobility vehicles such as the Mini Crosser with a maximum speed of up to 10 km/h count as cycles in traffic legislation and the traffic regulations for cycles must be followed at speeds in excess of 6 km/h.

At speeds of up to 6 km/h the Mini Crosser is regarded as a "pedestrian", and it is permissible to travel on the pavement at this speed. "Tortoise" speed corresponds to approx. 6 km/h. Lights must only be used during lighting-up time. Crash helmet and seat belts are not compulsory.

## Electromagnetic compatibility



If the Mini Crosser starts making involuntary movements or if the brakes are released, turn the Mini Crosser off as soon as it is safe to do so. In certain circumstances a Mini Crosser can set off shop alarms.

The Mini Crosser satisfies the requirements for the use of scooters in an environment with electromagnetic noise. There may, however, be rare situations in which electromagnetic noise can affect the Mini Crosser. Sources of such noise include radio and television stations, amateur radio transmitters and mobile phones.

If such equipment is being used close by, it is recommended that the Mini Crosser is switched off. Mobile phones should not be used while driving.

## When driving in traffic...

Be particularly aware of the following when driving in traffic:

- The Mini Crosser is a low vehicle and not always easy for other road users to see. Make quite sure that other road users have seen you before driving onto the highway.
- Keep an eye on traffic behind you. Keep well over to the side of the road when driving on busy roads.
- Turning right and left at crossroads. Be aware of cyclists and pedestrians. Follow the rules of the road for cyclists.
- How quickly things are happening. How long do the lights stay green? How quickly are cars approaching? etc.

## Specific driving situations

### Up kerbs

- Stop at right angles to the kerbstone about 5-10 cm away from it. Keep an eye on other road users.
- Lean forwards.
- Accelerate moderately so that your vehicle can overcome the obstacle. Do not stop halfway, but reduce speed once the vehicle is up.
- If the kerb is too high, do not try again, but find an alternative route.



### Down kerbs

- Lean back.
- If you are driving onto a road with traffic, keep an eye on other road users.
- Drive forwards and down the kerb at low speed. Make sure that your anti-tilt wheels (if fitted) do not catch on the edge.



### Up a ramp/hill

- There is a risk of tipping over backwards if the seat is pushed back when you start driving up a slope, e.g. a ramp.
- Pull the seat forward! Check that any ramp is stable.
- Lean forwards.
- Accelerate moderately so that your vehicle can overcome the obstacle. Do not stop halfway. Reduce speed once the vehicle is up. If you need to start on a hill, accelerate slowly so as not to tip over backwards.



### Down a ramp/hill

- Check that any ramp is stable.
- Lean back.
- Drive slowly down. Avoid stopping midway on short, steep slopes. On long hills, it is advisable to stop every now and again if you feel your speed is getting too high.



### **Across a slope**

- Lean into the slope.
  - Avoid sudden and sharp movements. This is particularly relevant when reversing.
  - Always drive at low speed.
- 

### **Anti-tilt wheels / stabilisers**

The Mini Crosser is a very stable vehicle. HOWEVER, in the case of incorrect weight distribution or careless driving there is a risk of tipping over. On the M Joy this risk is increased due to the more laid back position, and that there is no steering column to hold on to.

That is why the anti-tilt wheels are standard equipment on M Joy. (See picture below). In the front is mounted a weight block.



*Anti-tilt wheels*



## Batteries

The Mini Crosser uses sealed, maintenance-free GEL batteries (Exide 12 V/56 Ah or AGM 12 V/75 Ah).

They do not generate gas and do not have to be topped up with water.

---

## Charging



Please NOTE that the Mini Crosser can be equipped with several types of charger (ask your dealer for information on the various types).



**IMPORTANT!** Only ever use a charger designed for charging dry maintenance-free batteries. The max. charging current is 12 A.

If charging is to take place outdoors, an enclosed charger without a fan should be chosen.

The battery manufacturer recommends charging the batteries at temperatures between +10° and +30°C in order to achieve the charging times specified in “Technical data”.

Charging will take longer at temperatures below +20°C. This is because the battery finds it more difficult chemically to absorb the current.



Please NOTE that the capacity of the batteries will reduce over time and at low temperatures. In practice this means that vehicles with old batteries have a shorter driving distance than they had when the batteries were new.


Battery capacity at -10°C is about half the equivalent capacity at +20°C.

Avoid draining the batteries completely. We recommend charging the batteries in a heated room. If it is not possible to do this every day, it is a good idea to charge the Mini Crosser in a heated room (approx. 20°C) once a week.

New batteries do not reach full capacity until they have been charged and discharged about 20 times.



It is **IMPORTANT** to fit the batteries correctly. The battery terminals and snap locks are marked +/- . They must be fitted as shown in the sketch below. Make sure that the snap locks are properly closed. For the same reason there must not be any burrs on the terminals.



The Mini Crosser should be charged while not in use. The charger supplied from the factory switches over to trickle charging (very low power consumption) automatically when the batteries are fully charged. You can therefore leave the charger connected until the Mini Crosser is next used. The charger cannot overcharge the batteries.


The charger will flash until charging is complete. Then it will show a steady light.



The charger must not be placed on the seat during charging.



**IMPORTANT!** Never use charging devices other than those supplied from the factory without first contacting the dealer.



*NEUTRIK NC3MX charging plug*



*Connect the charging plug to the 3-pin socket on the joystick.*

## Battery disposal

Used batteries must be disposed of through your supplier or at a recycling centre.

Take care when handling leaky batteries, as the gel contains corrosive acid.

**INFO!**

New batteries can be purchased from Medema Production A/S.

---

## Storage

The scooter should be stored and charged under cover, preferably at temperatures above 0°C.

**Note!**



The charger must be kept dry, but should not be covered when in use.

In the case of long-term storage we recommend covering the Mini Crosser M Joy to protect it from dust, rain and sunlight.

---

## Cleaning

Clean the Mini Crosser M Joy with a damp cloth. If necessary, wipe dry with a chamois leather.

Wipe the covers dry and polish with car wax.



**IMPORTANT!** Never use a high-pressure cleaner or hose, as it could damage the Mini Crosser's electronics.

---

## Changing the wheels

If you get a puncture in one of the pneumatic tyres or if a tyre is so badly worn that it needs to be replaced, follow the instructions below.

Tyres and inner tubes can be purchased from the authorised dealer who supplied the Mini Crosser.

The Mini Crosser must be TURNED OFF before you start.



### Changing the wheels on the 4W model

5 mm Allen key

- Undo the five bolts.
- Take the wheel off.
- When the wheel is put back on, the spring washers must be refitted between the wheel rim and bolts.

The bolts must be done up tightly.



### REMEMBER!

- to let all the air out before taking the wheel apart to repair a puncture
- to secure the bolts with spring washers

## Fuses

### The Mini Crosser has the following protection systems:

- Thermal, which disconnects the power to the motor when overloaded. This thermal fuse is placed on the back of the scooter. The fuse protects against overloading of both steering and engine. When overloaded, it will switch off the engine - full speed and traction can be resumed after the motor has cooled down for 2-5 minutes.
- Mini Crosser M Joy also features a built-in fuse in the steering box. This will turn off, for example. hard driving in wet grass or very soft ground - wait 5 secs. Before you turn on the Mini Crosser again. This fuse will be activated before the thermal jump.



*Thermal fuse placement*

## Brakes

There are three braking systems on the Mini Crosser M Joy:

- Motor brake
  - Magnetic brake
  - Handbrake
- 
- Motor brake: adjusts vehicle speed also when going downhill.
  - Magnetic brake: automatic and comes on when accelerator is released. It also acts as a parking brake when the Mini Crosser is at a standstill. In an emergency, the Mini Crosser can be stopped instantaneously by turning of the joystick. Please note that this will cause very sharp braking. The rear wheels will lock!  
**Must not be used in the normal course of driving!**
  - Handbrake: intended as an emergency brake and parking brake. It must be operated with caution when driving in slippery conditions and downhill.

## Disengagement

- Turn the vehicle off on the joystick (I/O).
- Push the disengagement lever at the back down.  
The motor brake is now disengaged and the vehicle can be pushed or towed, but not driven.



### NOTE!

**The motor brake must never be disengaged on sloping terrain.**

**Once the motor brake has been disengaged, the Mini Crosser can only be braked with the hand brake.**



*Push the disengagement lever down to disengage (only handbrake works) and push up to reactivate the motor (normal driving now possible).*

*Disengagement lever*

---

## Transporting by motor vehicle

The Mini Crosser must always be restrained and the handbrake locked during transport in a motor vehicle or trailer.

Avoid lifting by the seat, covers, handlebars and armrests

If you just need to lift the scooter slightly, take hold of it between the rear lights and by the front bumper.

Secure it in the vehicle with belts attached to the two eyes at the front and two at the back. All the eyes are marked in yellow. See the section entitled “Securing to vehicle floor with belts”.


## Securing to vehicle floor with belts

Dahl Engineering belt set for securing in motor vehicles. Item No. C2-0242


ALWAYS use four belts at the back and two at the front.

The belts must always be attached to approved fittings in the vehicle and the four eyes welded to the scooter.

The belts MUST be attached within the angles shown in the picture for optimum security.









For transport in an estate car, the Mini Crosser must be secured with belts to the floor of the car


*Eye for restraints*

The Mini Crosser can be used as a seat during transport in a van or bus, provided it is securely anchored to the vehicle with approved four-point belts attached to the restraint points provided on the Mini Crosser.

The Mini Crosser's restraint points have been tested and approved in accordance with ISO 7176-19.

The user must also always be independently restrained in the actual motor vehicle in accordance with traffic legislation rules.

### Example



#### Restraining the passenger with a static 3-point seat belt:

Secure to the rearmost retractors.

The shoulder belt must rest against the collar bone and fall diagonally to the hip, where it is secured.

Tighten the belt by pulling on the loose strap. Undo it again by lifting the buckle. This is the same as on an aircraft.



Retractor with belt



Male and female parts of belt done up



Remember to turn the Mini Crosser off during transport. (I/O button).

However, if at all possible, we recommend that the user occupies one of the seats in the motor vehicle. All other things being equal, this is safer.

---

## Transporting by plane

If the Mini Crosser is to be transported by plane, the airlines require:

- the batteries to be flight-approved
- the air to be let out of the tyres
- the battery leads to be disconnected (not always, but frequently)

The seat and cover have to be removed to disconnect the battery leads.

A battery declaration for air travel can be found on the Mini Crosser website:

[http://www.minicrosser.dk/Download\\_brochurer.asp](http://www.minicrosser.dk/Download_brochurer.asp)

---

## Towing



If you should be unfortunate enough to break down, the Mini Crosser can be towed or pushed. The Mini Crosser must always be turned off and the motor brake disengaged during towing. See the section on Brakes.



If the Mini Crosser is to be towed, secure a rope to the tow fitting on the front - marked with a yellow "hook mark". Do not tow faster than 5 km/h. The scooter will generate electricity when it is towed, with the motor acting as a dynamo. If it is towed at more than 5 km/h, there is a risk of the motor generating enough electricity to damage the scooter and, in the worst case, cause a fire.

---

## Dimensioned drawing, Mini Crosser M Joy 4W



## Introduction to seats - Ergo Standard



### Lever for seat rotation

Pull the lever back to release the seat. The seat can then be turned 90° to either side with a stop every 45°.

The lever is spring-loaded and the seat locks automatically when the lever is released.

This lever is located on the right as standard, but can also be put on the left.



### Lever for seat forward/back

Pull the lever out/up to release the seat on the slide rail. The seat can then be moved forward or back as required.

When the lever is released, the seat automatically locks at the nearest stop.

There is a chock block on the front and back of the slide rail.



### Adjusting the angle of the armrest

The adjusting screw makes it possible to alter the angle of the armrest by about 15°.

The arm can be folded right up to make it easier to get in and out.

(Height adjustment of the armrest is an optional extra.)




### **Adjusting the angle of the backrest**

The angle of the backrest can be adjusted by about 30°. This can be done by undoing the adjusting screw and changing the position of the spacer. Do the screw up again after making the adjustment.




### **Optional extra:**

The seat can be fitted with a gas cylinder for adjusting the angle of the backrest.



### **Armrest forward / back:**

The armrest can be moved forward / back on the C-rail. This can be done by undoing the screws indicated. Do the screws up again after making the adjustment.



### **Backrest forward / back:**

The depth of the seat can be adjusted. This can be done by moving the backrest forward and back while the seat cushion is held in position. Undo the screws indicated on both sides. The backrest can

then be moved as far back as the edge of the rail and about 15 cm forward. Do the screws up after making the adjustment. The screws must be at least 2 cm from the end of the C-rail.

## Introduction to seats - Mando



### Adjusting the angle of the backrest:

The angle of the backrest can be adjusted by about 30° by raising this lever.




### Adjusting the angle of the armrest:

The angle of the armrest can be adjusted by turning this screw.



### Adjusting the height of the armrest:

The armrest can be raised/ lowered by undoing the two screws and moving the bracket one hole up or down.



### Headrest:

The headrest can be raised and lowered steplessly.

## Troubleshooting

The following is a list of various problems that, in our experience, may occur. The list gives possible causes and remedies.



Problem	Possible causes	Remedy
The Mini Crosser will not go.  The battery indicator is not lit.	The joystick has not been turned on (I/O).  The batteries are completely flat. The control fuse has blown. The main fuse have blown.	Turn on the joystick (I/O) and wait 5 sec. before activating the accelerator. Charge the batteries.  Change the fuse. Contact supplier.
The Mini Crosser will not go, but the battery indicator is lit.	The Mini Crosser has been overloaded.  The handbrake is on. There is a fault in the electronics. The batteries are flat. The charging plug has not been removed.	Wait approx. 1 min. before trying again. The vehicle must be turned off ( <a href="#">see section on Fuses</a> ). Release the handbrake. Contact supplier.  Contact supplier. Remove the charging plug.
The driving speed is too low.	The speed selector is on slow. The electronics are overloaded. There is too little air in the tyres.	Change to a faster speed.  Stop and wait a few seconds before starting. Pump the tyres up to the right pressure.
The driving distance per charge is too short.	There is a problem with the batteries. Low temperature.  There is a problem with the charger. There is too little air in the tyres. The charging method is wrong.	Charge the batteries and check that the green lamp on the charger lights up before driving off. Contact supplier.  Pump the tyres up to the right pressure. Read the section on <a href="#">Charging</a> in the User Manual.
The charging lamp on the charger does not light up when the charger is connected to the mains and the Mini Crosser.	No power to the switch. Fault in cable. Fault in charger.	Turn the switch on. Contact supplier. Read the operating instructions for the charger. Contact supplier.
The "ready" lamp on the charger does not light up even though the charger has been on for 10-12 hours.	There has been a power cut. The charger is doing a top-up charge. There is a problem with the batteries. There is a fault in the charging plug for the Mini Crosser M Joy.	Reconnect the charger and repeat the charging process. Check again half an hour later. Contact supplier.  Push the charging plug all the way in and repeat the charging process. Read the operating instructions for the charger.





Problem	Possible causes	Remedy
The "ready" lamp on the charger lights up even when partly discharged batteries are connected.	The fuse in the charger has blown. The switch in the charging plug is malfunctioning.	Contact supplier. Contact suppliers - read the operating instructions for the charger.
The charger lamp is showing an error.	The charging plug has not been inserted or there is a mains fault. The battery voltage is too low for charging to start.	Push the charger plug in or contact the supplier. Read the operating instructions for the charger - or contact the supplier.

## Troubleshooting the Dynamic

If there is an electronic fault, a number of lamps in the battery indicator will light up. The following table shows what they mean.

 <p><i>Status light</i></p>	 <p><i>Status light</i></p>
<b>Troubleshooting the G90T and G90A</b>	
Status light flashes with 1 flash per cycle	Switch off the scooter and restart. If this does not work, contact a competent service engineer. (HMC or therapist).
Status light flashes with 2 flashes per cycle and the scooter runs at half speed.	The scooter may be fitted with a speed limiter that halves the speed when the seat is raised. This is not a fault.
Status light flashes with 3 flashes per cycle.	Fault, left motor. Contact a competent service engineer. (HMC or therapist).
Status light flashes with 4 flashes per cycle	Fault, right motor. Contact a competent service engineer. (HMC or therapist).
Status light flashes with 5 flashes per cycle.	Fault, left motor brake. Contact a competent service engineer. (HMC or therapist).
Status light flashes with 6 flashes per cycle	Fault, right motor brake. Contact a competent service engineer. (HMC or therapist).
Status light flashes with 7 flashes per cycle.	The batteries need charging. Or they have a defect.
Status light flashes with 8 flashes per cycle.	The battery voltage is too high.

	
<p>Status light flashes with 9 flashes per cycle.</p>	<p>Electronic fault. Contact a competent service engineer. (HMC or therapist).</p>
<p>Status light flashes with 10 flashes per cycle.</p>	<p>Electronic fault. Contact a competent service engineer. (HMC or therapist).</p>
<p>Status light flashes with 11 flashes per cycle.</p>	<p>The motor has exceeded the maximum voltage, or has come close to exceeding it for longer than the parameter value. Contact a competent service engineer. (HMC or therapist).</p>
<p>Status light flashes with 12 flashes per cycle.</p>	<p>Wrong combination of modules. Contact a competent service engineer. (HMC or therapist).</p>
<p>The scooter runs at half speed or not at all.</p>	<p>If the cause is not a raised seat, this may be a serious fault. Contact a competent service engineer as soon as possible. (HMC or therapist).</p>



The Dynamic electronic system is programmable so that the driving characteristics can be optimised for individual users. Programming must only be carried out by people trained in Dynamic's control systems. Changing parameters incorrectly can produce driving characteristics that are dangerous for the user.



**Note**

Mobile phones should not be used very close to the scooter as they generate electromagnetic fields that may affect the electronic system. The electronic system of the scooter itself also generates some electromagnetic fields, which on rare occasions may affect equipment like the alarm systems in shops.




Dynamic's control systems have been tested and meet the requirements of ISO7176/14 and EN12184


## Programming



NOTE! For safety reasons, modifications may only be made by trained personnel, such as Mini Crosser service engineers and consultants or authorised service personnel at mobility centres, etc.

## Technical data

Technical data	M Joy
<b>General information:</b>	
Transport height without seat	109 cm
Total length: Total length incl. anti-tilt wheel: Total width:	140 cm 148,5 cm 68 cm
Total weight incl. batteries and Ergo Std seat.	152 kg
Weight without seat (Ergo Standard 45 cm).	131 kg
Dynamic stability in all directions.	10° - 17%
Kerb climbing	11 cm
Max. speed	10 km/h
Braking distance 10 km/h	2.0 metres
Turning radius	141 cm
Max user weight - standard	150 kg
<b>Wheels</b>	
M-127-3-165 13x5.00-6" ext. Ø325 mm Wheel Norway black - T, E, Nordic, MaxX, MaxX HD and M-model	4,1 bar 
M-127-3-165P 13x5.00-6" ext. Ø325 mm Wheel Norway black w/spikes - T, E, Nordic, MaxX, MaxX HD and M-model	4,1 bar 
M-127-3-065 13x3.00-8" ext. Ø340 mm Wheel black - T and M-model	3,5 bar 

Technical data	M Joy
T-127-3-068 13x3.00-8" ext. Ø340 mm Wheel black puncture free - T and M-model	PUR 
Classification, ISO	Class C
<b>Seats:</b>	
Effective seat width: Ergo Standard Spinalus Mando HD (user weight max 250 kg) Ergo child / junior	40,45,50,55,60 cm 40, 45, 50 cm 48 cm 60, 65, 70 cm 35, 37 cm
Effective seat depth: Ergo, Spinalus and HD seats Mando Ergo child / junior	32-52 cm 48 cm 20-30 and 25-40 cm
Seat height: Ergo standrad / HD Spinalus Mando Ergo child / junior	48 cm 56 cm 53 cm 36-41 cm
Height: front edge of seat to ground Ergo Standard with slide rail Ergo Standard without slide rail Mando seat with slide rail Mando seat without slide rail	64 - 74 cm 61 - 71 cm 64 - 74 cm 61 - 61 cm
Height: footplate to front edge of seat Ergo Standard with slide rail Ergo Standard without slide rail Mando seat with slide rail Mando seat without slide rail	45 - 55 cm 41 - 51 cm 45 - 55 cm 41 - 51 cm
Angle of backrest Ergo standard / HD (manual) Ergo standard / HD (EL) Spinalus Mando Ergo child / junior	Approx. -5° to +20° Approx. -5° to +30° Approx. -5° to +20° Approx. +5° to +30° Approx. -5° to +20°
Backrest hight	48 cm
<b>Batteries:</b>	
56 Ah batteries	Standard
Maximum driving distance with new batteries at +20° on flat, firm surface	35 km
75/80 Ah batteries	Option
Maximum driving distance with new batteries at +20° on flat, firm surface	55 km

<b>Technical data</b>	<b>M Joy</b>
110/115 Ah batteries	Option
Maximum driving distance with new batteries at +20° on flat, firm surface	75 km
Optimum battery capacity is reached after approx. 20 charges / discharges.	
Battery type: Standard: Option: Option:	2 x 12 V / 56 Ah 2 x 12 V / 75 Ah 2 x 12 V / 110 Ah
Max. battery dimensions in cm.	26.5 x 17 x 22 28.4 x 26.7 x 23
Battery weight, 2 batteries 56 Ah 75/80 Ah 110/115 Ah	43 kg 50 kg 80 kg
Energy consumption in kWh, when charging from "empty" 56 Ah 75/80 Ah 110/115 Ah	Approx. 1.5 Approx. 1.5 Approx. 3.0
Charging device, 24 V DC 56 Ah 75 Ah 110 Ah	6 - 10 A 6 - 10 A 10 - 12 A
Approx. charging time at 20° C	8 hours
<b>Lighting:</b>	
Bulb, headlight	24 V - 21 W
Diode, rear light	24 V - 2 W
Diode, indicator	24 V - 2 W
Standard colour	Grey - metallic

## International addresses

### Australia / New Zealand

#### Pride Mobility Products

Tel: +61 3 9706-4611  
Fax: +61 3 9706-4622  
Email: prideaustralia@pride-mobility.com.au  
Internet: www.pridemobility.com/international/Australia/australia.asp

### Belgien

#### Mobile Belgium

Tel: +32 (0)50 317919  
Fax: +32 (0)50 311026  
Email: info@mobilescooter.be  
Internet: www.mobilescooter.be

### Denmark

#### Head Office/Jutland Medema Danmark A/S

Tel: +45 7010 1755  
Fax: +45 9718 8835  
Email: info@medema.com  
Internet: www.medema.dk

### Denmark

#### Zealand

**Medema Danmark A/S**  
Tel: +45 7010 1755  
Fax: +45 5664 3388  
Email: info@medema.com  
Internet: www.medema.dk

### Denmark

#### Funen

**Medema Danmark A/S**  
Tel: +45 7010 1755  
Email: info@medema.com  
Internet: www.medema.dk

### Finland

#### Respecta OY

Tel: Sales and service +358 (0)207649748  
Email: myynti@respecta.fi  
Fax: +358 (0)207649735  
Tel: Maintenance - parts and service: +358 (0)207649747  
Email: huolto@respecta.fi

### France

#### ICARE

Tel: +33 (0)1 60 78 94 84  
Fax: +33 (0)1 69 13 09 78  
Email: icarejmb@wanadoo.fr  
Internet: http://www.icare-jmb.com

### Germany

#### Baden Württemberg Bavaria

#### Draisin GmbH

Tel: (+49) 0 78 41 66 77 0  
Fax: (+49) 0 78 41 66 77 88  
Email: info@draisin.com  
Internet: www.draisin.com/

### Germany

#### Baden Württemberg

#### Elofahr Palm

Tel: +49 07225 918 831  
Fax: +49 07225 988 070  
Mobile: +49 0171 314 7726  
Email: uwepalm@web.de

### Germany

#### Berlin, Mecklenburg-West Pomerania, Saxony-Anhalt, Thuringia, Brandenburg and Saxony

#### Conze Elektro-Mobile

Tel: +49 030 6789 7969  
Fax: +49 030 6789 7968  
Email: info@elektro-mobile.de  
Internet: www.elektro-mobile.de

### Germany

#### Hesse, Rhineland-Palatinate, Saarland and Lower Franconia

#### Beromobil GmbH

#### Elektromobile Treppenlifte

Tel: +49 06123 - 60 57 22  
Fax: +49 06123 - 60 57 23  
Email: info@beromobil.de  
Internet: www.beromobil.de

### Germany

#### Schleswig-Holstein, Lower Saxony and North Rhine- Westphalia

#### Wiggers GmbH - Reha / Rollstuhl Werkstatt

Tel: (+49) 0441 36 111 302  
Fax: (+49) 0441 36 111 309  
Email: wiggers.gmbh@ewetel.net  
Internet: www.minicrosser.info; www.ot-wiggers.de

### Greece

#### AMEASHOP

Tel.: 210 7777545  
Fax.: 210 7777835  
Mobile: 6970021890  
Internet: www.ameashop.gr  
Email: faidros@ameashop.gr

## Holland

**Richard van Seenus Nederland B.V.**  
Tel: +31 (0) 36 53 20 450  
Fax: +31 (0) 36 53 21 308  
Email: info@vanseenus.nl  
Internet: www.vanseenus.nl

## Island

**Icepharma HF**  
Tel: +354 540 8000  
Fax: +354 540 8001  
Email: icepharma@icepharma.is  
Internet: www.icepharma.is

## Italy

**Comfort Online s.r.l. - Società Unipersonale**  
Tel.: +39 041 5461105  
Fax: +39 041 5462098  
Email: info@comfortonline.it  
Internet: www.comfortonline.it

## Norway

**Medema gruppen AS**  
Tel: +47 6706 4900  
Fax: +47 6706 4990  
Email: firmapost@medema.no  
Internet: www.medema.no

## Poland

**J.A.S.**  
Tel.: +48 667 432 282, +48 605 674 017  
Fax: +48 33 810 54 58  
Email: biuro@jas-bielsko.pl  
Internet: http://jas-bielsko.pl/

## Spain

**Assessorament I Mobilitat. S.L.**  
Tel: +34 93 666 55 23  
Fax: +34 93 666 88 85  
Email: info@assessoramentimobilitat.com  
Internet: http://www.assessoramentimobilitat.com

## Sweden

**Minicrosser AB**  
Tel: (+46) 303 245 200  
Fax: (+46) 303 245 228  
Email: info@minicrosser.se  
Internet: www.minicrosser.se

## Switzerland

**Degonda Rehab SA**  
Tel: +41 033 345 83 00  
Fax: +41 033 345 83 01  
Email: info@turbo-twist.ch  
Internet: www.turbo-twist.ch

## UK

**Mini Crosser Ltd**  
Tel: Sales +(44) 0800 043 0029  
Tel: Service & Warranty +(44) 0800 169 9588  
Email: sales@mini-crosser.co.uk  
Internet: www.mini-crosser.co.uk

## UK

**Active Mobility**  
Tel.: +(44) 0135 77 44 22  
Internet: www.active-mobility.co.uk

## USA / Canada

**Pride Mobility Products Corp.**  
Tel: +1 800 800 8586  
Fax: +1 800 800 1636  
Internet: www.pridemobility.com



**medemagroup**